

GoLift LED Functions

GoLift
compact ceiling lift



Charging	Battery	LED	Buzzer	Additional Notes
YES	Full	Green If re-connect to charger LED goes Orange for 5 seconds then back to Green	N/A	LED stays on 1 minute after use then shuts off
NO	Full	Green	N/A	LED stays on 1 minute then shuts off
YES	Not Full	Orange	N/A	LED stays on 1 minute then shuts off
NO	Low Battery	Blinks Orange	N/A	LED stays on 1 minute after use then shuts off. Monitor voltage to check for battery level. One full lift at max load must be able to be completed.
NO	Dead	Red	Single beep every 3 seconds	LED stays on 1 minute after use then shuts off. Monitor voltage to check for battery level. Lift in up direction not allowed but must be able to go in the down direction.
NO	Fault Temperature Error	Flashing Red	Beep pattern single beep every 5 seconds	LED stays on for 1 minute
Emergency Lowering		Red	Beep	Red LED if main board is working No LED if main board has malfunctioned

Latch Error: Red flashing once every 5 second and buzzer on solid

During Sleep Mode (when LED turns off - after 1 minute) the unit should draw no power (1ma). This is required to ensure long shelf life. When any buttons are pressed the unit has to wake up instantly.

Maintenance Alarm

- At 1001 lifts, the LED will flash Green
- To reset, press UP button, then the DOWN button, then UP button, then DOWN button and then hold both buttons for 10 seconds
- If successful, the LED turns Red, Green and Orange. The buzzer beeps three times
- Please note, that a reset can be performed at any time, not just after the 1000th lift.

For further assistance with your GoLift Hoist, contact:

sales@wealdenrehab.com

01634 813388



Wealden Rehab
Care Equipment Excellence

wealdenrehab.com

0845 658 8411