

RAZ ATP Shower Chair Assessment Form



1 Therapist and Users Details

Assessment Date:

Therapist Details

Name:

Organisation:

Email:

Telephone:

User Details

Name/Ref:

Address:

.....

.....

Email:

User Weight:

Notes:

.....

.....

i Information
If you require assistance in completing this form, please call our sales office on 01634 813388,
or use our online chat facility at www.wealdenrehab.com/chat



2 Core Stock Chair

RAZ ATP Standard frame:

Frame Code N0001P

plus



N676

Standard Armrests



N650P

Visco Foam Seat
457 x 280 - 406mm



N721

MFX Footrests



N731

Standard Headrest



N2581

Commode pan

Ordering details to include:

- Seat Height
- Armrest Height
- Backrest Angle



Information

Seat measurements can be configured on **pages 6 & 7**.

Set Up:

On delivery the prescriber is required to adjust:

- Tension Adjustable Backrest
- Seat Depth
- Footplate Height




2 Chair Dimensions

1

Seat Width

Measure the outside of the users pelvis




- ☐ 450mm
- ☐ 400mm (achieved with hip guides)
- ☐ 350mm (achieved with hip guides)
- ☐ 300mm (achieved with hip guides)
- ☐ Other

2

Backrest Width


Measure shoulder to shoulder



- ☐ 355mm
- ☐ 400mm
- ☐ 450mm
- ☐ Other

3

Backrest Height




- ☐ 380mm
- ☐ 460mm
- ☐ Other

4

Seat Depth

Measure user's back to back of knee



- ☐ 280mm
- ☐ 300mm
- ☐ 330mm
- ☐ 360mm
- ☐ 380mm
- ☐ 410mm
- ☐ Other



2 Chair Dimensions

RAZ ATP Seat

5

Standard Seat



- N650P
- Visco Foam Seat for RAZ ATP
- Bridged Front
- 457 x 280 - 406 mm Adjustable
- Tapered Oblong Aperture with 64mm thickness

Other Seat Options



Visco Foam Open Front Seat



Visco Foam Side Access Seat



Visco Foam Rear Open Seat

Custom seats available on request




2 Chair Dimensions

Seat Height

6

When choosing the seat height consider clearance required for the toilet, side transfer or standing transfer heights.



Toilet Bowl Height: (seat up)

.....

Chair Seat Height:

.....


Toilet bowl height (mm)	380	400	420	440	460	480	500	520
ATP Standard Seat Height	530	530	560	560	580	610	640	640



2 Chair Dimensions

Lining up your chair with the toilet bowl

Rearmost chair point to rear of seat aperture



ATP = 230mm

AP = 180mm

SP = 280mm

AT = 290mm


Using an automatic WC?

Please contact Wealden Rehab to discuss options.

7

Footplate Height

Measure user's lower leg length



Seat to Footrest dimensions:

200-300mm, 300-400mm and 400-500mm

i Information

If the standard footplates are not suitable - see the accessories list of footplates on **page 8**.



3 Accessories

Footplate Accessories*



N722 – MFX Footrests for RAZ Rehab Shower Chair, pair, minimal forward extension, adjustable 300-410mm. For users with a lower leg length between 300-410mm.



N721 – MFX Footrests for RAZ Rehab Shower Chair, pair, minimal forward extension, adjustable 200-300mm. For users with a lower leg length between 200-300mm.



N720 – Standard Footplates, pair, Adjustable 41-52cm, with folding plate.



N724 – Angle Adjustable Footplates for RAZ Rehab Shower Chair. Footplate can be angled to accommodate users with footdrop or if the ankle is in a fixed position.



N2628 – Neoprene Footplate Covers, for RAZ Standard Footplates. To give a softer texture under the users foot.



N2632 – Calf strap to give the legs support when the chair is tilted and to prevent feet from slipping backwards.

*Additional charge



3 Accessories

Postural Support Accessories*



N2676P – Long PU pad for use with Leg Adductor Brackets – for RAZ-AT (Pair). This is provided to give side support on the chair. If a user is abducting, this will help their legs on the seat cushion.



N690P – Swing-away Lateral Supports for RAZ Shower Chair, pair, complete with adjustable lateral pad and neoprene support wrap. These give thoracic support to the patient's torso, helping to keep their posture upright. They can be adjusted in height/depth and angle.



N692 – Lateral Support Extension Plates for RAZ Rehab Shower Chair. This gives 25mm additional adjustment per side (50mm in total for a pair). These are used to reduce the width between the lateral support or the adductor pads.

Height from seat to top of lateral:

.....

Distance between laterals:

.....

Headrest Support Accessories*



N731 – Standard headrest.



C5960 – Head support cushion. This head cushion is flexible so can be curved around the users head. Comes with a stem to give height and depth adjustment.



C5961 – Head support cushion. Very solid head support for clients who cannot support their own head. This support is not flexible.



C5963 – Head support cushion. Most common headrest for users needing lateral support for their head. This is flexible so can wrap around the sides of the head to help support it in an upright position.



N3906 – RAZ Headrest Clamp with 75mm Side-To-Side Offset Adjustment – specify this if the support required is not in the centreline of the chair.

*Additional charge



3 Accessories

Chair Accessories*

Backrest



N2620 – Reclining AdjustaBack for RAZ Shower Chair, 500mm height, with adjustable flip-up armrests. Reclines to 30°. NB Chassis Tilt will be limited to 25°.

Arms



N676 – Standard Armrests Flip-up for easy access and transfers (come as part of Adjustaback)



N676 – Armrest Lock Kit for RAZ Rehab Shower Chair, pair. This kit stops the armrest moving up and therefore can be used for the patient to pull on and re-position. The lock releases to allow the arm to lift freely.



N681 – These arm troughs help to keep the arm on the pad without falling off. There is angle and depth adjustment to the pads to suits the position of the users arms.



N679 – Wide Soft Armpad for RAZ Rehab Shower Chair, 350mm x 80mm – Pair. This armrest offer a wider area of support for the patient's arms.



N683 – Pivoting Arm Mount – Pair. These are used in conjunction with product codes N681 and N679, in order to allow the armpad to rotate through 90° and move into the upright position without conflicting with the backrest.

Self Propelling Wheels



N3157 – Self propel kit for AT. A kit to change an attendant tilt chair to a self propelled tilt chair. This kit fits to any standard RAZ-AT and includes anti-tips and brake levers. Wheels with textured push rim.

Commode Pan



N2581 – Commode Pan with handle and lid. Autoclavable.

Chest Harness



Flex-Contol harnesses offer a greater range of posture
5.71032 – Mens
5.71034 – Womens

Lap Belts



3001 - Shower Hip Belt, Single Adjustable



6.70042 -4-point Hip Belt, Single Adjustable

*Additional charge




4 Sizing Guide


Dimension Chart – ATP




Max Patient Weight
115kg




A Overall Width
507mm




B Width between handles
360 / 410 / 460mm




C Height from top of seat to top of backrest vinyl
280/460mm




D Height from floor to top of push handles*
990mm




E Top of seat from floor
530 / 560 / 580
610 / 640mm




F Underside of chassis side cross bar
180mm




G Height of arm above seat
200 – 260mm



H Depth of chassis
450mm




I Overall depth without footrests
890mm




J Overall depth with footrests
980mm




K Top of seat from standard footplates
410 – 520mm




L Rearmost chair point to rear of seat aperture
230mm



M Front of chassis to front of footplates
250mm



N overall depth with footplates, in full tilt
1130mm



O maximum clearance for WC cistern*
825mm



5