



# Saros Modular Chair

## Assessment Form

# 1 Therapist and Users Details

Assessment Date:

## Therapist Details

Name: .....

Organisation: .....

Email: .....

Telephone: .....

## User Details

Name/Ref: .....

Address: .....

Email: .....

User Weight: .....

Pressure Care: .....

Posture: .....

Mobility: .....

Diagnosis: .....

Prognosis: .....

Transfer Method: .....

Current Seating: .....

Notes: .....

## Information

If you require assistance in completing this form, please call our sales office on 01634 813388, or use our online chat facility at [www.wealdenrehab.com/chat](http://www.wealdenrehab.com/chat)

## 2 Saros Chair Specifications

The Saros chair provides a comfortable chair on castors that is adaptable for different patient sizes and users that require pressure relief or additional support to sit with a good posture.

- Easily accessible push handles for carers
- Multi-adjustable pillow backrest as standard
- Integrated pressure care options
- Multi tilt range
- Multi back angle adjustments
- Arm height adjustment
- Adjustable seat width
- Multi seat ramp angle
- Webbed base seating system
- Footplate
- Four braked castors



### Fabric Care

Wipe clean routinely with a microfibre cloth, lightly dampened with cold water only. Then dry with a clean absorbent cloth.

Soak up spillages as soon as possible with an absorbent dry cloth to prevent staining. Clean with a microfibre cloth, lightly dampened with cold water.

Note dyes and pigments from indigo jeans, ball point and felt tip pens may permanently stain.

If soiling has accumulated use a mild soap (non-alkaline). Take care to remove all soap residues. Then dry with an absorbent cloth.

Remove body fluids and extreme soiling immediately (within 15 minutes) with cold water. Apply your approved infection control protocol, a 10,000ppm sodium chlorite bleach solution, Haz-Tabs or Chlor-Clean. Remove product residues and then dry with a clean absorbent cloth.

Note steam cleaning disinfects but the application of heat may fix staining substances and may result in permanent marks.

Do not apply solvent based cleaners or any other chemicals to the fabric.


Do not dry clean, tumble dry or fully immerse the fabric in water.

## 3 Chair Configuration

### Seat Width

1

Measure user's  
hip to hip width




Choose closest seat width

- ☐ 405mm/16"
- ☐ 457mm/18"
- ☐ 508mm/20"

### Seat Depth

2

Measure user's base  
of spine to back of knee




Choose closest seat depth

- ☐ 432mm/17"
- ☐ 457mm/18"
- ☐ 482mm/19"
- ☐ 508mm/20"
- ☐ 533mm/21"
- ☐ 559mm/22"

### Seat Height

3

Measure user's  
seat height



Choose seat height between  
405-432mm / 16"-22"

4

Choose backrest angle



90°

STANDARD



100°



110°

5

Choose seat ramp angle



Flat

STANDARD



5°



10°

### Information

If the user's dimensions are outside of these parameters, please contact Wealden Rehab.

## 4 Pressure & Posture Management

### Cushion choice

6

#### Visco-elastic cushion



Visco-elastic cushion – High risk pressure relieving cushion, to assist in prevention of pressure sores.



STANDARD

#### Integrated cushion



Cut-out allows integration of an alternating air cushion for users with existing pressure sores.



OPTIONAL

### Backrest choice

7

#### Butterfly backrest



To provide extremely high levels of comfort, designed to accommodate patients with spinal curvature or in support where there may be muscle weakness.



STANDARD

#### Curved thoracic wedges



Pair of curved wedges to provide additional thoracic support, often of benefit for users with neurological conditions.



OPTIONAL

#### Thoracic backrest support



Thoracic backrest support – For users with kyphosis, offering a head cushion that protrudes further forward for support.



OPTIONAL

### Headrest choice

8

#### Profiled horseshoe headrest



Fits over the shoulders for users with reduced head control.



OPTIONAL

# 5 Accessories

## Additional Support

9

### 90° Footplate



90° fixed angle padded footplate.



STANDARD

### Angle Adjustable Footplate



Gives support to feet and reduces risk of foot drop. Vapour permeable fabric.



OPTIONAL

### Pelvic Strap



Safe, comfortable positioning aid.



OPTIONAL

### Four Point Pelvic Strap



Fixes to chair at four points for additional pelvic control for safe, comfortable positioning aid.



OPTIONAL

## Actuation

10

### Manual Operation



Fully manually operated chair.



STANDARD

### Electric Operation



Battery powered operation of tilt-in-space and legrest elevation allows the user to adjust the chair independently.




OPTIONAL


## 6 Dimensions Chart




Max Patient Weight  
150kg / 23.5st




Overall Width  
65-75cm / 26-30"




Overall Height  
124cm / 49"




Overall Depth  
105-127cm / 42-50"




Seat Width  
40-50cm / 16-20"




Seat Depth  
40-55cm / 16-22"



Arm Height  
7", 8", 9"



Back Height  
76cm / 30"



Back Angle  
90, 80, 70, 60°

7